# ZERO NANOTECH

## Ion Microscopy, Machining, and Elemental Analysis with the Cesium Low Temperature Ion Source (LoTIS)



Adam V Steele adam@zeroK.com

Visit us at <u>www.zeroK.com</u>

August 2022

# Technology and Applications









- Low Temperature Ion Source
  - Laser-cooling + Photoionization
- Heavy ion nanomachining
- Small spot sizes
- Excellent resolution at low energy (~2 nm resolution at 1 pA, 16 kV)
- 1 pA 10 nA

#### FIB:ZERO

- LoTIS + FIB
- Comparable to standard Ga<sup>+</sup> FIB, with 2x higher resolution at low beam currents
- Compatible with normal peripherals, gas chemistries etc..



#### SIMS:ZERO

- FIB:ZERO with SIMS
  - Analysis of secondary ions in a mass spectrometer
- Best for elemental-compositional analysis
- Collab. with Luxembourg Institute of Science and Technology (LIST)

### FIB:ZERO





# Milling Homogeneity: 150 nm Au on Si $\rightarrow$ Cs<sup>+</sup> LoTIS proves even touchdown

Milled with Ga<sup>+</sup> LMIS



### Milled with Cs<sup>+</sup> LoTIS







- milled rectangle 'almost through' the Au layer
- milling time Ga and Cs almost the same

# Milling Accuracy: 110 nm Au on Si $\rightarrow$ LoTIS provides clean mill boxes with sharp corners

### Milled with Ga<sup>+</sup> LMIS



# TECHNISCHE UNIVERSITÄT KAISERSLAUTERN



Milled with Cs<sup>+</sup> LoTIS

- squares with 1, 0.6, 0.4, 0.2, 0.1 and 0.05 μm length
- milled through the Au layer
- milling time Ga and Cs almost the same

### FIB:ZERO Cross Section Example



Done with 200 pA beam, 30 pA 'cleanup' afterwards



### FIB:ZERO Cross Section Example

Oversaturated images to show lack of curtaining



ZERØK

# Depth of Focus Comparison

(Results on slides that follow)

### "Wood Pile" Nanostructures

- Heights: 40 μm, 80 μm, 120 μm
- In the following slides we acquire an image containing both the top and bottom of such the 120 µm (tallest) structure
- We can compare the depth of focus of various beams by comparing the 'blurriness' of the top of the structure

A better depth of focus aids in the milling and imaging of 'deep' or 'tall' structures.

			50			80.46 µm (cs)				121 <u>.5 µm (c)</u>		Nano Strue
*	HV	curr	dwell	det	mode	WD	tilt	mag 🎛	HFW	200 μm —		

### FEI: SEM image



TECHNISCHE UNIVERSITÄT

### Depth of Focus Comparison →LoTIS depth of focus substantially better than Ga

### Ga<sup>+</sup> LMIS (30 kV)

Cs<sup>+</sup> LoTIS (10 kV)



"Wood Pile" Height 120  $\mu m$ 

ZERØK

TECHNISCHE UNIVERSITÄT KAISERSLAUTERN

# Invasiveness Comparisons (SRIM Calculations)

- Comparison of three scenarios where spot size might be 'good enough'
- Cs has significantly reduced straggle and implant depth



# FIB:ZERO Milling Rates



Milling rate of 10 kV Cs<sup>+</sup> FIB:ZERO about 15% lower than 30 kV Ga<sup>+</sup> for Si

Cs<sup>+</sup> LoTIS milling rates 90% higher than Ne<sup>+</sup> (and **much** higher than He<sup>+</sup>)

Ne 10 kV	Ga 30 kV	Cs 10 kV
1.00-1.38 at/ion	2.20-2.40 at/ion	1.90-2.15 at/ion

# SUMMARY FIB:ZERO

#### ... is a 'nanomachining' tool

... has industry-leading performance at low beam currents and low energy

... is compatible with gas precursors for deposition or etch just like other FIBs

Data pictured right, implant depth and milling fidelity, summarize the story best

#### Invasiveness Comparisons (SRIM Calculations)

8

ZERØK

Comparison of three scenarios where spot size might be 'good enough'
Cs has significantly reduced straggle and implant depth



Confidential- DO NOT DISCLOSE



#### 12

# Pain Points of Existing Elemental Analysis Techniques and a New Solution



### EDX/EELS

- Long sample-prep times
- 3D analysis infeasible
- Low-Z elements challenging

### Site-Specific SIMS

- Resolution limited to ~50 nm with high yield (CAMECA NanoSIMS), or
- Can get a high resolution FIB (Ga, He, Ne) with a time-of-flight SIMS analyzer. But low secondary ion yields from these beams usually results in poor lateral resolution. Additionally, time-of-flight analyzers necessitate **long** acquisition times.

These points are addressable by **SIMS:ZERO** 

- Few-nanometer resolution (slide 21)
- High secondary ion yield (slides 23,24)
- Integrated sample-prep and analysis capability (slides 25-31)

These slides are some of the first public data presented from SIMS:ZERO

# SIMS:ZERO

Focal Plane

(4X)

### Instrument Overview

Cs+ FIB:ZERO (zeroK) and SIMS spectrometer (LIST: Luxembourg Institute of Science and Technology) on a 600 series FIB (FEI)

Secondary Ion Beam Axis (+ or -) - FIB online 6/2020

- SIMS online 5/2021





### Primary Ion Species Matters



Differing Sputter Rates  $\rightarrow$  Analysis Time

Differing interaction Volumes  $\rightarrow$  Resolution

Differing Yields → Sensitivity Floor, SNR





# FIB / SIMS Combination

### Sample Prep, Nanofabrication / Analysis, Process Control







#### **LoTIS** capabilities

- 2-16 keV Cs+ beam
- Up to 5nA beam current
- Spotsize <2nm at low current
- Good spotsizes even at low beam energy

#### FIB Mode (SIMS Extraction Optics Retracted)

- Milling
- Sample Preparation (eg Sectioning, Polishing)
- Nanofabrication
- Gas-assisted processes (eg Platinum Deposition)
- Tilt stage

#### SIMS Mode (SIMS Extraction Optics Inserted)

- Highest spatial resolution SIMS imaging
  - $\sigma = 6 \text{ nm}$  demonstrated
- Mass resolution  $M/\Delta M = 400$
- Mass range up to 300 amu
- High secondary ion throughput (~40% simulated)
- 4-Channel Detector Standard (Continuous Focal Plane Detector available)

ZERØK

### SIMS:ZERO Resolution Tungsten Carbide



- SIMS:ZERO can provide higher resolution SIMS scans than any other instrument
- SIMS resolution is a function of abundance, yield, and spot size
- SIMS:ZERO has a focused ion beam with <3 nm spot size, and since it's Cs<sup>+</sup> we achieve high yields for many materials
- In samples with high abundances, resolution at near the physical limits of SIMS can be achieved (see right)



Multi\_WC\_2105121624015\_CH1.TIF



### $\sigma = 6.1 nm$ (!)

272s acquisition time.

Working Distance = 51.6mm

Sample	WC (184 amu)
FOV (um)	2.97um
I (pA)	2.5
U (kV)	16

### SIMS Analysis Example CIGS Cu(In,Ga)Se<sub>2</sub> – Rb doped



### <u>Summary</u>

- CIGS is a solar cell absorber material
  - Rubidium doping increases conversion efficiency
- SIMS spectra clearly show all CIGS elements:
  - Cu, In, Ga, Rb in Positive Mode
  - Se in Negative Mode
- Secondary ion imaging channels show distribution of elements in sample, eg Rb dopants concentrated in grain boundaries
- Secondary electron images provide complementary information at high resolution
- Section view technique provides superior SIMS data





Werner, et al. <u>Scientific</u> <u>Reports</u> volume 10, 7530 (2020)

### CIGS Cu(In,Ga)Se<sub>2</sub> – Rb doped Mass Spectra – Positive Ions





Start (mT): 30.000000 Stop (mT): 700.000000 Delta (mT): 0.100000 Sampling rate (ms) : 250.000000 Waiting time (s) : 0.250000 Period of beam acq:0 Pos CH1: 100.000675 Pos CH2 : 200.000362 Pos CH3: 299.999717 Pos CH4: 390.000000 Date 06/23/2021 CIGS Sample Aperture Slit 100um 43 FOV (um) 10 I (pA) U (kV) 16

### CIGS Cu(In,Ga)Se<sub>2</sub> – Rb doped Mass Spectra – Negative Ions





BScan Neg CIGS 2106241202174.csv 105 C CN Se CH3 104 ក្ល 10<sup>3</sup> 102 10<sup>1</sup> 20 40 60 80 100 120 140 160 180 200 220 240 260 0 amu

Start (mT): 30.000000 Stop (mT): 700.000000 Delta (mT): 0.100000 Sampling rate (ms) : 250.000000 Waiting time (s) : 0.250000 Period of beam acq:0 Pos CH1: 100.000675 Pos CH2: 200.000362 Pos CH3: 299.999717 Pos CH4: 390.00000 Date 06/24/2021 Sample CIGS Aperture Slit 100um FOV (um) 43 I (pA) 10 U (kV) 16

# SIMS-Compatible Section View 45° Angle Cut - Example





View with Sample Normal to Beam; Ready for SIMS on Section



For many samples, working with a section view is a sensible choice

- 1. Reveal sub-surface structure
- 2. Obtain depth profile data without accumulated topography from uneven sputtering
- Polish rough samples to isolate elemental from topographical contrast
- 4. Build 3D tomographic reconstructions through serial sectioning/polishing

In SIMS:ZERO, sample must be normal to ion beam in SIMS Mode, so section face is cut at 45° to sample surface

### CIGS Cu(In,Ga)Se<sub>2</sub> – Rb doped Serial Sectioning / Imaging / Polishing Work-Flow



SE Images



SIMS section, prepared with low surface topography, reveals layer structure (glass, moly, CIGS, Window/Buffer Stack) After SIMS Imaging, section face develops topography which obscures elemental contrast / distribution information Section face after cleanup mill. Ready for SIMS on next layer

Cs+, 16keV, 10pA, 51.6mm WD

### CIGS Cu(In,Ga)Se<sub>2</sub> – Rb doped Section View – Positive Ions



#### SE Image – Pre-SIMS



- Rb confined to grain boundaries
- Grains are smaller near the interfaces

Cs+, 16keV, 3.5pA, 51.6mm WD CIGS\_Pos\_2107161606287.csv CIGS\_Pos\_2107161613425.csv

### CIGS Cu(In,Ga)Se<sub>2</sub> – Rb doped Section View – Negative Ions – Post 2<sup>nd</sup> Polish





SE Image – Post Polish Low topography restored





Signal band in CIGS layer near moly may be sulfur, commonly used in CIGS fabrication process; inclusions near surface



Se is more uniformly distributed in CIGS layer; droplets at moly interface, a few inclusion near surface

> Cs+, 16keV, 10pA, 51.6mm WD CIGS\_Neg\_2107201513310.csv

### CIGS Cu(In,Ga)Se<sub>2</sub> – Rb doped Section View – Positive Ions – Post 3<sup>rd</sup> Polish





Na – Soda Lime Glass

Ga concentration gradient ↑



In concentration gradient  $\downarrow$ 

Cs+, 16keV, 10pA, 51.6mm WD CIGS\_Pos\_2107201626359.csv

### CIGS Cu(In,Ga)Se<sub>2</sub> – Rb doped Section View – Positive Ions





#### Apparent width of Rubidium signal between grains



Cs+, 16keV, 3.5pA, 51.6mm WD CIGS\_Pos\_2107151409368.csv

## Diatoms **LIST** (Silica — shelled algae) 7.5 um FoV





0







# Location of TiO nanoparticle within in huge, fixed cell <sup>12</sup>C<sup>14</sup>N











LIST

# Application Example: SIMS:ZERO as EDX Alternative



EDX elemental analysis is capable of few-nm resolution and can image the majority of elements well, but sensitivity is limited to a few tenths of a percent and sample prep is time consuming

Historically, SIMS has offered excellent (ppm) sensitivity but limited lateral resolution

Now, SIMS:ZERO enables creation of elemental maps with both few-nm resolution and excellent sensitivity without lamella preparation

These capabilities also make possible the creation of 3D elemental maps

**Existing** Workflow - Thin Sample EDX



Only one shot : analysis limited to a single depth

#### **Optimized** Workflow - SIMS:ZERO

Load Sample SIMS Analysis **SIMS** Analysis **SIMS** Analysis Locate Prepare SIMS Analysis on Deeper on Deeper on Deeper ROI **Cross Section** SIMS:ZERO Plane Plane Plane

### Continuous Detector A SIMS:ZERO Option

- SIMS signals for a given element are split into many lines (e.g. Ti, TiO etc..).
- In discrete-detector systems this leads to a loss of information and lower SNR.
- With continuous detector technology we can sample the entire mass spectrum at once.
- Now we can collect the entire spectrum as in TOF systems, but without painfully long acquisition times.

(Top) Photo of SIMS spectrometer (at LIST) and the continuous focal plane detector mounted to a vacuum flange.(Bottom) A 480mm micro-channel plate that spans the focal plane of the spectrometer.

TiO<sub>2</sub> NPs on InP wafer

1000

100

2000

2500 3000

3500

4000 Channel

4500

sdo

113In

5000 5500





LIST



## SUMMARY SIMS:ZERO

... has all the capabilities of FIB:ZERO

... adds high-resolution, high-sensitivity, high speed elemental analysis

... enables new modes of operation

Slides at right featuring new modes of operation and high resolution elemental mapping summarize the story best

### CIGS Cu(In,Ga)Se<sub>2</sub> – Rb doped

Section View – Positive Ions

SE Image – Pre-SIMS







ZEROK

- Rb confined to grain boundaries
- Grains are smaller near the interfaces



### Application Example: SIMS:ZERO as EDX Alternative

EDX elemental analysis is capable of few-nm resolution and can image the majority of elements well, but sensitivity is limited to a few tenths of a percent and sample prep is time consuming

FOV9.5

listorically, SIMS has offered excellent (ppm) sensitivity but limited lateral resolution

Now, SIMS:ZERO enables creation of elemental maps with both few-nm resolution and excellent sensitivity without amella preparation

These capabilities also make possible the creation of 3D elemental maps

Existing Workflow - Thin Sample EDX



Only one shot : analysis limited to a single depth

Optimized Workflow - SIMS:ZERO

AVS67, October 2021

## Spot Sizes

### Selected Beam Energies and Currents

Results given as a  $\sigma$ 

- $d_{50} = 2.2 \sigma$
- $d_{35-65} = \frac{\sigma}{1.3}$ ,
- $d_{16-84} = 2\sigma$

Spot sizes are about 2x smaller than Helios (Ga<sup>+</sup>) at <10 pA, and at lower energy

Bit worse spot sizes than Zeiss He/Ne, but better machining performance in many cases

Methodology in [1]

[1] A V Steele *et al* 2017 *Nano Futures* **1** 015005 Open access link: <u>https://iopscience.iop.org/article/10.1088/2399-1984/aa6a48</u>

16	kV	8 kV			
Current (pA)	Spot Size (1- $\sigma$ nm)		Current (pA)	Spot Size (1- $\sigma$ nm)	
1.5	<2.0		1.5	2.3	
3.0	2.3		3	2.5	
10	4.0		10	4.7	
30	7.5		30	7.6	
100	23		100	55	
300	57		300	150	
1000	175		1000	244	
5500	580		2600	510	

### Low-Energy

Energy (kV)	Current (pA)	Spot Size (1- $\sigma$ nm)
5	3.5	15
2	3.5	44



# → FIB:ZERO can provide very narrow metal lines

#### Pt lines with 100nm width as deposited on Si.





#### Tilted view to show height of Pt lines above Si substrate.

